


*NanoScope™ Signal Access Module (SAM)—  
Description & Use*

---



# *NanoScope™ Signal Access Module (SAM)— Description & Use*

## *210.1 Description*

The NanoScope™ Signal Access Module (SAM)<sup>1</sup> enables SPM users to make direct connections to their microscope’s electronics. For most applications, the SAM avoids disassembling the microscope, and is provided with BNC connectors which interface with most standard laboratory instruments (see Figure 210-1). Twenty-five separate SPM signals are available. Toggle switches enable the user to switch each line separately between a normal, uninterrupted configuration and an external input signal. Output BNC connectors permit monitoring of both conditioned and uninterrupted signals for use in experiments.



**IMPORTANT NOTE:** While providing maximum access to SPM electronics, the SAM also makes the unit more vulnerable to abuse. Users are cautioned that Digital Instruments assumes no responsibility for repair costs to the SPM or peripheral equipment associated with misuse of the SAM.

Signal connectors are divided into three groups:

- Low Voltage Inputs to NanoScope Controller (data signals)
- Low Voltage Outputs from NanoScope Controller (control signals)
- High Voltage Outputs from NanoScope Controller (piezo drive signals)



<b>Rev.</b>	<b>Date</b>	<b>Section(s) Affected</b>	<b>Ref. DCR</b>	<b>Approval</b>
Rev. C	29MAY97	Section 210.2	171	<i>Lulan</i>
Rev. B	05JUN96	New formatting	N/A	<i>BA</i>
Rev. A	26OCT95	Released	N/A	<i>BB</i>

1. The Signal Access Module functions as a “breakout box.”

In addition, the SAM also features BNC connectors for  $\pm 15$  VDC, +5 VDC and an RS232 9-pin plug. At the rear of the module are located two D-37 connectors for attachment to the microscope and NanoScope controller

Normally, the SAM is connected between the control box and microscope. If an Extender™ Electronics Module is utilized, the SAM is connected between the Extender and the controller.



**Figure 210-1 Signal Access Module (SAM)**

---

The underside of the SAM features silk-screened information vital to safe, successful use of the module and should be reviewed by all users. **DO NOT EXCEED MAXIMUM RATED VOLTAGES/CURRENTS OR THE MICROSCOPE AND CONTROLLER MAY BE DAMAGED!** This information is provided in Appendix A of this support note.

## Safety Precautions

### Danger—High Voltage!:

The bottom row of the SAM features seven output lines with maximum voltages of 220 VAC and currents up to 70 ma. The unit is shipped with all high voltage BNC connectors covered by red insulating caps. **DO NOT REMOVE INSULATING CAPS UNTIL BNC CONNECTORS ARE TO BE USED! RECOVER CONNECTORS NO LONGER IN USE.** Do not insert foreign objects (screwdrivers, etc.) into high voltage BNC connectors. Accidental contact with the high voltage contacts may result in serious injury.



Users are advised to disconnect the entire unit before any connections are made to high voltage lines. Be especially cautious whenever handling fluids on or around the microscope. For those imaging samples in fluid, it is recommended that fluids be kept in spill-proof squirt bottles and used sparingly. Wipe away all spilled fluids immediately and avoid spilling fluids into the microscope's imaging head.

### Protect External Equipment


A 1.0 amp fuse is wired into the 5 VDC line; all other lines are unfused. If protection of external equipment is required, it is recommended that fuses or optoisolators be utilized.

When making multiple connections between the SAM and external equipment, it is recommended that all connection cables be labeled clearly to prevent confusion. External equipment which is being modified to operate with the microscope should be modified *before* connection to the SAM.

Test all new connections and setups with a voltmeter or oscilloscope *before* connecting to the microscope to verify minimum and maximum output levels. Output levels should be compared with values listed on the underside of the box (see Appendix A) to verify that they are within safe operating limits.

## Lift Up on Toggle Switches

Toggle switches used with the SAM are designed to avoid accidental toggling—they must be pulled upward to toggle, then released (see Figure 210-2). This “locks” the toggle switch in position. **DO NOT FORCE TOGGLE SWITCHES!**



**Figure 210-2** Lift switches when toggling.

## Do Not Remove Base

Do not remove the base of the SAM. There are no user-serviceable parts inside and a risk of high voltage shock is present. If failure of the SAM is suspected, return the entire unit to Digital Instruments for repair or replacement.

## 210.2 Installing the SAM

1. Remove the SAM from its packing box. The packing box should contain the following items:


- SAM
- Connector Cable, D-37 type, 4-foot.
- Insulator Covers (7), red. Installed on all high voltage outputs.
- *Support Note No. 210* (this document).

2. Place the SAM in a convenient, level location between the microscope and the NanoScope controller.

3. Verify that all toggle switches on the SAM are toggled to “Output” (toggled right), then power down the SPM and NanoScope controller.

**NOTE:** Switches should be toggled to “Output” (right) during normal, uninterrupted operation.


4. If you have an Extender module the labels on the SAM will be correct if you place the SAM between the NanoScope controller and the Extender box (refer to figure below).



**Figure 210-3 “Normal” Configuration**

---

If you place the SAM between the Extender and the SPM, some of the lines will be different as noted.



**Figure 210-4 Alternate Configuration**

<u>Lines</u>	<u>Alternate<sup>1</sup> Signals</u>
AuxA	Lateral deflection
AuxD	STM preamp output
Ana1	Piezo drive signal
Ana2	Tip/Sample voltage (as selected by switches on outside of Extender Module)
D1	D1 (if switch inside Extender Module is set for Dimension) RMS (if switch is set for MultiMode)

5. Repower the SPM and NanoScope controller.

1. Note: Please refer to Appendix A for the “normal” configuration signals.



### 210.3 Making Connections to the SAM

**IMPORTANT:** Under normal conditions, external devices should not interfere with imaging electronics. However, to ensure the integrity of image data, verify that any device which is used to insert a signal into the SAM has an output impedance of  $\leq 1\text{ K}\Omega$ . Devices receiving signals from the SAM (e.g., oscilloscope) should have input impedances of  $\geq 1\text{ M}\Omega$ . If it is necessary to work beyond these impedance limits, contact Digital Instruments for instructions.

#### 210.3.1 Monitoring Signals from the SAM

The SAM may be used to monitor any signal from the NanoScope controller, or any signal measured or produced by the microscope. To monitor a signal from the SAM, do the following:

1. Verify that each monitoring device to be connected has an impedance of  $1\text{ M}\Omega$  or greater.
2. Make certain that all switches are toggled to the right position, except when an external signal is being inserted into the SAM (see "Putting Signals into the SAM" below).
3. Connect all monitoring devices to the appropriate outputs using shielded BNC connectors and cables. Signals are labeled below each BNC pair. Refer also to Drawing No. ACC-BOB-8154.




**CAUTION!** High voltages are directly accessible from the SAM (bottom row). *Do not touch high voltage contacts!* Before connecting, verify that all monitoring devices to be used are compatible with the outputs by checking Appendix A or the underside of the SAM. Improperly connecting NanoScope high voltage signals can damage both the NanoScope and external devices.



### Example A. Monitoring a Feedback Signal (In 0)

The below diagram shows how to make a simple connection for monitoring feedback signals using an oscilloscope; the connection is made at the In0 connector on the upper-left corner of the SAM. The switch is toggled to the right position.




#### 210.3.2 Putting Signals Into the SAM


There are two reasons for connecting external signals to the SAM: 1) to insert data signals into the NanoScope III controller; 2) to provide an external control signal or piezo drive signal to the microscope. See Examples B—F below.

1. Verify that input devices to be connected have output impedances of less than 1 K $\Omega$ , and voltage/current outputs that do not exceed ratings provided on the underside of the SAM and in Appendix A.
2. Connect the external signal to the appropriate “Input” BNC connector. (Signals are identified below each pair of connectors and on Drawing No. ACC-BON-8154.)
3. Toggle each selected signal’s line switch to the left, “Input” position. **Note:** toggling any switch to the left will interrupt the standard electrical connection.
4. Use inputs while staying within voltage/current limits rated on the underside of the SAM and in Appendix A. As inputs are added, closely monitor image data for any abrupt changes. If abrupt voltage changes are noted and/or if the image fails suddenly, disconnect the input line immediately—the input line may be faulty.


**Example B. Connecting an External Data Signal to the SAM**




**Example C. Connecting an External Piezo Drive Signal**




**Example D. Connecting an External Control Signal**



**Example E. Modifying a Control Signal with an External Signal (Modulating the Bias Voltage)**



### Example F. Modifying a Data Signal (E.g., Lock-in Amplifier on Feedback Signal In 0)



#### 210.4 Cautions

Drawing No. ACC-BOB-8154 (attached) is a schematic diagram of the SAM, and illustrates its internal connections. Users may refer to this drawing when connecting external equipment. Users are advised to test all new setups using a voltmeter or oscilloscope before connecting to the SAM. Input levels should remain within limits listed on the underside of the SAM and in Appendix A.

The signal access module provides unlimited access to the NanoScope's input and output signals. This ensures maximum flexibility for customization; however, it also provides ample opportunity to damage the controller, microscope and external devices through misuse. For example, connecting a high-voltage output to a low-voltage digital output will surely destroy one or both of these output lines. To protect against such mishaps, please carefully follow the following guidelines:

1. Always check and recheck signal labels before making external connections.
2. Carefully check voltage, current and impedance limits before connecting any external device.
3. Never connect two outputs directly together.

## 210.5 Appendix A

**Signal Access Module Connector Ratings**

SIGNAL	MAX. VOLTAGE	MAX. CURRENT	SIGNAL TYPE
In0	± 10 VDC		Analog Input
AuxA	± 10 VDC		Analog Input
AuxB	± 10 VDC		Analog Input
AuxC	± 10 VDC		Analog Input
AuxD	± 10 VDC		Analog Input
Zmod	± 10 VDC		Analog Input
A	± 5 VDC	500 mA	OCL Output
B	± 5 VDC	500 mA	OCL Output
C	± 5 VDC	500 mA	OCL Output
D	± 5 VDC	500 mA	OCL Output
D0	± 5 VDC		TTL Output
D1	± 5 VDC		TTL Output
Ana1	± 10 VDC	10 mA	Analog Output
Ana2	± 12 VDC	10 mA	Analog Output
Ana3	± 10 VDC	10 mA	Analog Output
Ana4	± 10 VDC	10 mA	Analog Output
Bias	± 10 VDC	10 mA	Analog Output
LV Z	± 12 VDC	10 mA	Analog Output
X	± 220 VDC	70 mA	High Voltage Output
$\bar{X}$	± 220 VDC	70 mA	High Voltage Output
Y	± 220 VDC	70 mA	High Voltage Output
$\bar{Y}$	± 220 VDC	70 mA	High Voltage Output
Z	± 220 VDC	70 mA	High Voltage Output
$\bar{Z}$	± 220 VDC	70 mA	High Voltage Output
Ana2 (HV)		70 mA	High Voltage Output
+ 15 VDC		500 mA	Power Supply
-15 VDC		500 mA	Power Supply
+ 5VDC		500 mA	Power Supply

RS232 Connector: Pin 2 = RXD; Pin 3 = TXD; Pin 5 = GND.

Appendix A, Support Note 210, Rev. C Microscope Product-Mode Signal Matrix


LINES		MultiMode™						Dimension™-series SPMs & BioScope™					STM			Contact AFM							
		TappingMode	Contact AFM	ECAFM	STM	ECSTM	STS	Force Modulation	TappingMode	Contact AFM	STM	STS	Force Modulation	Standard	ECSTM	STS	Standard	LFM	Tip View STM	STS	ECAFM		
Analog Inputs to NanoScope® ± 10VDC	In 0	RMS amplitude	Vertical deflection <sup>a</sup> minus setpoint voltage		Voltage as a function of tunneling current				RMS amplitude	Vertical deflection minus setpoint voltage	Voltage as a function of tunneling current			Voltage as a function of tunneling current			Vertical deflection <sup>a</sup> minus offset voltage	Vertical deflection <sup>b</sup> minus offset voltage	Voltage as a function of tunneling current		Vertical deflection minus offset voltage		
	Aux A	Vertical deflection <sup>a</sup>	Lateral deflection (LFM) <sup>a1</sup>	cell current	Set by user	cell current			Vertical deflection <sup>a</sup>	Lateral deflection (LFM) <sup>a1</sup>				Set by user			Set by user	Lateral deflection (LFM) <sup>b1</sup>	Set by user		Cell current		
	Aux B	Set by user			Set by user			See note <sup>c</sup>													Not used		
	Aux C	Set by user			Set by user			Photodiode array sum (A+B+C+D)															
	Aux D	Set by user. See note <sup>e</sup>		cell potential	See note <sup>e</sup>	cell potential			Set by user. See note <sup>d</sup>													Cell potential	
Analog outputs from NanoScope® ± 10 VDC ± 220 VDC	Bias	Setpoint voltage			Sample bias				Setpoint voltage		Sample bias			Sample bias			Setpoint voltage		Sample bias		Setpoint voltage		
	Ana.1	Cantilever drive amp.	Set by user	Working electrode potential	Set by user	Tunneling tip potential			Cantilever Drive amp.	Set by user				Set by user			Set by user				Working electrode potential		
	Ana.2 LV	Configured at factory to control 1X & 8X output attenuation of final stage, but may be reconfigured using jumpers in MultiMode base. <sup>e</sup>												Set by user			Set by user						
	Ana.2 HV	22 X Ana.2 LV. See note <sup>e</sup>												Set by user									
Digital Outputs from NanoScope®	D0	See note <sup>f</sup>	line sync	cell on/off	line sync	cell on/off			See note <sup>f</sup>	line sync				line sync	Cell on/off		line sync				Cell on/off		
	D1		frame sync	galv/pot	frame sync	galv/pot				frame sync				frame sync	galv/pot		frame sync				galv/pot		
Motor Lines from NanoScope	A	write strobe/stepper motor	Stepper motor						write strobe/stepper motor	Stepper motor						Stepper motor				Stepper motor			
	B	serial in/stepper motor						serial in/stepper motor															
	C	frequency load/stepper motor						frequency load/stepper motor															
	D	reset synth./stepper motor						reset synth./stepper motor															

## Appendix A, Support Note 210, Rev. C Microscope Product-Mode Signal Matrix

a. Vertical deflection (see photodiode array-laser beam orientation below):  $\frac{(A + B) - (C + D)}{A + B + C + D}$

a'. Lateral (horizontal) deflection (see photodiode array-laser beam orientation below):  $\frac{(A + C) - (B + D)}{A + B + C + D}$


Photodiode orientation (MultiMode™, Dimension™-series SPMs, standard contact AFM, MFM, EFM):



b. Vertical deflection— Contact AFM Lateral Force Microscope (LFM):  $\frac{A - B}{A + B}$

b'. Lateral (horizontal) deflection— Contact AFM Lateral Force Microscope (LFM) friction measurements:  $\frac{C - D}{C + D}$

Photodiode orientation (LFM head only):




c. Some auxillary lines in Dimension-series SPMs are used for checking limit switches. Previous to engagement, the system alternates between TappingMode and contact AFM so that the LFM signal may be used to update the photodiode display.

d. Aux D connects to a jumper on the backplane that is normally jumpered to STM. Aux D may also serve as a user input.

## Appendix A, Support Note 210, Rev. C Microscope Product-Mode Signal Matrix

e. MultiMode™ microscopes feature an internal interconnect board (P/N MMD-INT-1293) having jumpers. These jumpers are accessible through the opening in the bottom cover of the MultiMode base and may be reconfigured as shown here:



**NOTE:** Ana.2 HV (high voltage) is permanently set to be at 22 times (22X) the Ana.2 LV (low voltage)voltage.

f. The function of motor lines to the MultiMode base is controlled by D0 and D1 in the following ways:

<u>D0</u>	<u>D1</u>	Function
0	0	Motor lines control the engage motor.
0	1	Motor lines control the frequency synthesizer.
1	X	Motor and frequency synthesizer data is latched and motor lines from the NanoScope are ignored.

**NOTE:** When motor lines control the frequency synthesizer, functions listed on chart above for TappingMode motor lines (A = write strobe; B = serial in; C = frequency load; D = reset synthesizer) are in effect.



Appendix A, Support Note 210, Rev. C Microscope Product-Mode Signal Matrix

LINES		Interferometric AFM Detection					
		Manual Stage		StandAlone™		Large Sample Stage (LSS)	
		Contact AFM	TappingMode	Contact AFM	TappingMode	Contact AFM	TappingMode
Analog Inputs to NanoScope® ± 10VDC	In 0	detector					
	Aux A	Set by user					
	Aux B						
	Aux C						
	Aux D						
Analog outputs from NanoScope® ± 10 VDC  ± 220 VDC	Bias	Ground					
	Ana.1	Setpoint voltage					
	Ana.2 LV	Set by user					
	Ana.2 HV	Set by user					
Digital Outputs from NanoScope®	D0	line sync					
	D1	frame sync					
Motor Lines from NanoScope	A			Unused			
	B						
	C						
	D						

# Appendix A, Support Note 210, Rev. C Microscope Product-Mode Signal Matrix

## Magnetic Force Microscopy (MFM) and Electric Force Microscopy (EFM) Settings

LINES		MultiMode™				Dimension™-series SPMs & BioScope™				Interferometric AFM Detection												
		Without Extender™ Electronics Module		With Extender™ Electronics Module		Without Extender™ Electronics Module		With Extender™ Electronics Module		Manual Stage		StandAlone™				Large Sample Stage (LSS)						
		MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	
Analog Inputs to NanoScope® ± 10VDC	In 0	RMS amplitude				RMS amplitude																
	Aux A	Vertical deflection				Vertical deflection <sup>a</sup>																
	Aux B	Set by user.																				
	Aux C					Photodiode array sum (A+B+C+D)																
	Aux D	Set by user. See note <sup>e</sup>																				
Analog outputs from NanoScope® ± 10VDC  ± 220VDC	Bias	Setpoint voltage				Setpoint voltage																
	Ana.1	Cantilever drive amplitude				Cantilever drive amplitude																
	Ana.2 LV	Configured at factory to control 1X & 8X output attenuation of final Multi-Mode stage, but may be reconfigured using jumpers in MultiMode base <sup>e</sup> .																				
	Ana.2 HV	22 X Ana 2LV voltage.																				
Digital Outputs from NanoScope®	D0	See note <sup>f</sup>																				
	D1	See note <sup>f</sup>																				
Motor Lines from NanoScope	A	Write strobe/Stepper motor				write strobe																
	B	Serial in/Stepper motor				serial in																
	C	Frequency load/Stepper motor				frequency load																
	D	Reset synthesizer/Stepper motor				reset synth.																

Appendix A, Support Note 210, Rev. C Microscope Product-Mode Signal Matrix


LINES		MultiMode™						Dimension™-series SPMs & BioScope™					STM			Contact AFM						
		TappingMode	Contact AFM	ECAFM	STM	ECSTM	STS	Force Modulation	TappingMode	Contact AFM	STM	STS	Force Modulation	Standard	ECSTM	STS	Standard	LFM	Tip View STM	STS	ECAFM	
Analog Inputs to NanoScope® ± 10VDC	In 0	RMS amplitude	Vertical deflection <sup>a</sup> minus setpoint voltage		Voltage as a function of tunneling current				RMS amplitude	Vertical deflection minus setpoint voltage	Voltage as a function of tunneling current			Voltage as a function of tunneling current			Vertical deflection <sup>a</sup> minus offset voltage	Vertical deflection <sup>b</sup> minus offset voltage	Voltage as a function of tunneling current		Vertical deflection minus offset voltage	
	Aux A	Vertical deflection <sup>a</sup>	Lateral deflection (LFM) <sup>a1</sup>	cell current	Set by user	cell current			Vertical deflection <sup>a</sup>	Lateral deflection (LFM) <sup>a1</sup>				Set by user			Set by user	Lateral deflection (LFM) <sup>b1</sup>	Set by user		Cell current	
	Aux B	Set by user			Set by user			See note <sup>c</sup>													Not used	
	Aux C	Set by user			Set by user			Photodiode array sum (A+B+C+D)														
	Aux D	Set by user. See note <sup>e</sup>		cell potential	See note <sup>e</sup>	cell potential			Set by user. See note <sup>d</sup>													Cell potential
Analog outputs from NanoScope® ± 10 VDC ± 220 VDC	Bias	Setpoint voltage			Sample bias				Setpoint voltage		Sample bias			Sample bias			Setpoint voltage		Sample bias		Setpoint voltage	
	Ana.1	Cantilever drive amp.	Set by user	Working electrode potential	Set by user	Tunneling tip potential			Cantilever Drive amp.	Set by user				Set by user			Set by user				Working electrode potential	
	Ana.2 LV	Configured at factory to control 1X & 8X output attenuation of final stage, but may be reconfigured using jumpers in MultiMode base. <sup>e</sup>												Set by user			Set by user					
	Ana.2 HV	22 X Ana.2 LV. See note <sup>e</sup>												Set by user								
Digital Outputs from NanoScope®	D0	See note <sup>f</sup>	line sync	cell on/off	line sync	cell on/off			See note <sup>f</sup>	line sync				line sync	Cell on/off		line sync				Cell on/off	
	D1		frame sync	galv/pot	frame sync	galv/pot				frame sync				frame sync	galv/pot		frame sync				galv/pot	
Motor Lines from NanoScope	A	write strobe/stepper motor	Stepper motor						write strobe/stepper motor	Stepper motor						Stepper motor				Stepper motor		
	B	serial in/stepper motor						serial in/stepper motor														
	C	frequency load/stepper motor						frequency load/stepper motor														
	D	reset synth./stepper motor						reset synth./stepper motor														

## Appendix A, Support Note 210, Rev. C Microscope Product-Mode Signal Matrix

a. Vertical deflection (see photodiode array-laser beam orientation below):  $\frac{(A + B) - (C + D)}{A + B + C + D}$

a'. Lateral (horizontal) deflection (see photodiode array-laser beam orientation below):  $\frac{(A + C) - (B + D)}{A + B + C + D}$


Photodiode orientation (MultiMode™, Dimension™-series SPMs, standard contact AFM, MFM, EFM):



b. Vertical deflection— Contact AFM Lateral Force Microscope (LFM):  $\frac{A - B}{A + B}$

b'. Lateral (horizontal) deflection— Contact AFM Lateral Force Microscope (LFM) friction measurements:  $\frac{C - D}{C + D}$

Photodiode orientation (LFM head only):




c. Some auxillary lines in Dimension-series SPMs are used for checking limit switches. Previous to engagement, the system alternates between TappingMode and contact AFM so that the LFM signal may be used to update the photodiode display.

d. Aux D connects to a jumper on the backplane that is normally jumpered to STM. Aux D may also serve as a user input.

## Appendix A, Support Note 210, Rev. C Microscope Product-Mode Signal Matrix

e. MultiMode™ microscopes feature an internal interconnect board (P/N MMD-INT-1293) having jumpers. These jumpers are accessible through the opening in the bottom cover of the MultiMode base and may be reconfigured as shown here:



**NOTE:** Ana.2 HV (high voltage) is permanently set to be at 22 times (22X) the Ana.2 LV (low voltage)voltage.

f. The function of motor lines to the MultiMode base is controlled by D0 and D1 in the following ways:

<u>D0</u>	<u>D1</u>	Function
0	0	Motor lines control the engage motor.
0	1	Motor lines control the frequency synthesizer.
1	X	Motor and frequency synthesizer data is latched and motor lines from the NanoScope are ignored.

**NOTE:** When motor lines control the frequency synthesizer, functions listed on chart above for TappingMode motor lines (A = write strobe; B = serial in; C = frequency load; D = reset synthesizer) are in effect.

Appendix A, Support Note 210, Rev. C Microscope Product-Mode Signal Matrix

LINES		Interferometric AFM Detection					
		Manual Stage		StandAlone™		Large Sample Stage (LSS)	
		Contact AFM	TappingMode	Contact AFM	TappingMode	Contact AFM	TappingMode
Analog Inputs to NanoScope® ± 10VDC	In 0	detector					
	Aux A	Set by user					
	Aux B						
	Aux C						
	Aux D						
Analog outputs from NanoScope® ± 10 VDC  ± 220 VDC	Bias	Ground					
	Ana.1	Setpoint voltage					
	Ana.2 LV	Set by user					
	Ana.2 HV	Set by user					
Digital Outputs from NanoScope®	D0	line sync					
	D1	frame sync					
Motor Lines from NanoScope	A			Unused			
	B						
	C						
	D						

# Appendix A, Support Note 210, Rev. C Microscope Product-Mode Signal Matrix

## Magnetic Force Microscopy (MFM) and Electric Force Microscopy (EFM) Settings

LINES		MultiMode™				Dimension™-series SPMs & BioScope™				Interferometric AFM Detection												
		Without Extender™ Electronics Module		With Extender™ Electronics Module		Without Extender™ Electronics Module		With Extender™ Electronics Module		Manual Stage		StandAlone™				Large Sample Stage (LSS)						
		MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	MFM	EFM	
Analog Inputs to NanoScope® ± 10VDC	In 0	RMS amplitude				RMS amplitude																
	Aux A	Vertical deflection				Vertical deflection <sup>a</sup>																
	Aux B	Set by user.																				
	Aux C					Photodiode array sum (A+B+C+D)																
	Aux D	Set by user. See note <sup>e</sup>																				
Analog outputs from NanoScope® ± 10VDC  ± 220VDC	Bias	Setpoint voltage				Setpoint voltage																
	Ana.1	Cantilever drive amplitude				Cantilever drive amplitude																
	Ana.2 LV	Configured at factory to control 1X & 8X output attenuation of final Multi-Mode stage, but may be reconfigured using jumpers in MultiMode base <sup>e</sup> .																				
	Ana.2 HV	22 X Ana 2LV voltage.																				
Digital Outputs from NanoScope®	D0	See note <sup>f</sup>																				
	D1	See note <sup>f</sup>																				
Motor Lines from NanoScope	A	Write strobe/Stepper motor				write strobe																
	B	Serial in/Stepper motor				serial in																
	C	Frequency load/Stepper motor				frequency load																
	D	Reset synthesizer/Stepper motor				reset synth.																